


3DActive POSTURE CORRECTOR



USER GUIDE

THANK YOU FOR CHOOSING THE 3DACTIVE POSTURE CORRECTOR

This guide is available for download at: www.3dactive.com/guides

3DACTIVE.COM

CONTENTS

SAFETY WARNING	1
CARE INSTRUCTIONS	2
ADJUSTING YOUR POSTURE CORRECTOR	3
HELP	5

SAFETY WARNING

Please seek professional help and advice from your Doctor for further assistance when using your Posture Corrector.

- If you feel pain whenever you wear your Posture Corrector, stop using it and seek professional advice on how to best wear it/fit it to your posture.
- Make sure your Posture Corrector is not too tight as this can worsen your posture or even cause injury.



Start with 20 minutes a day and gradually increase wearing time. Make sure your Posture Corrector is not too tight to ensure comfortable wearing.

CARE INSTRUCTIONS

Please wash the Posture Corrector and make sure it's completely dry before storing it for long periods of time. Do not tumble dry the Posture Corrector or apply bleach.



Hand wash



30°C/80°F



Do not bleach




Do not tumble dry

ADJUSTING YOUR POSTURE CORRECTOR



Put your Posture Corrector on like you would a backpack or jacket.




Hold the ends of each velcro strap and pull them to tighten.



Bring both straps forward, under your arms.



Secure each velcro strap tightly over the shoulder area.



Your Posture Corrector should feel snug and comfortable.

NEED HELP?

HELP | QUESTIONS | WARRANTY

We stand by our product with excellent support.



SUPPORT CENTER

www.3dactive.com/support



EMAIL US

support@3dactive.com



CALL US

US +1 (855) 789 9777 - Mon-Fri 9am-5pm (Pacific Time)

UK +44 808 178 4009 - Mon-Fri 4pm-12am (GMT)

3DActive

3DActive Posture Corrector

Goldhat Limited | Dalton House, 60 Windsor Avenue, London, SW19 2RR, UK